



**Lantheus
Medical Imaging**

331 Treble Cove Road
North Billerica, MA 01862

800.362.2668
www.lantheus.com

November 2018

Dear Valued Customer,

This letter is to inform you that Lantheus Medical Imaging will begin to supply the Hospira Eluant Charge Vial, item number H-10, (10 mL Saline), compatible with our TechneLite®, Technetium Tc 99m Generator, effective early December 2018. We are making this change to ensure we are best positioned to provide you with a reliable supply that consistently meets your nuclear medicine requirements.

The Hospira Eluant Charge Vial, Hospira 0.9% Sodium Chloride Injection, USP (Preservative Free) 10mL vial – 25 vials/pack with NDC number 0409-4888-10, will continue to be shipped according to your current shipment schedule. Your Lantheus customer service representative will fulfill your current customer 6 vials/pack configuration (A17) requirements with the Hospira 25 vials/pack configuration that best matches your current eluant demand. As stated, the Hospira 10 mL saline vial is compatible with our TechneLite generator.

We recognize this change in supply requires a change in your practice. We sincerely apologize for any inconvenience this new supplier transition may cause and we appreciate your patience and support during our transition. Please know that we value your business, and are working diligently to transition to our new supplier as smoothly and as safely as possible.

Please note: If you requested and have received our Saline Charge Vial Alignment Guide, Lantheus recommends that you do not use the Saline Charge Vial Alignment Guide with the 10 mL Hospira saline vials. The Hospira 10 mL vial does fit into the Guide but, due to dimensional differences, the vial becomes flush with the top of the Guide when the saline vial is installed onto the saline vial charge well needle assembly.

Lantheus remains dedicated to providing a reliable supply of quality products, ensuring patient safety and meeting the nuclear medicine needs of our customers and the patients you serve. If you have additional questions, please contact your Customer Service Representative at (800) 299-3431.

Sincerely,

Dorothy Barr
Vice President, Global Supply Chain
Lantheus Medical Imaging, Inc.