
Turning Insights into Hope

May 11, 2010

Dear Radiopharmacist:

Spectrum Pharmaceuticals, Inc. received FDA approval on April 8, 2010 to use an additional supplier of Yttrium-90 (Y90) for the Zevalin® therapeutic regimen.

The additional supplier is Eckert & Ziegler Nuclitec, Inc. of Hopkinton, MA using their Auriga Medical™ branded Y90 produced at their cGMP manufacturing facility in Braunschweig, Germany.





Nuclitec

Beginning June 1, 2010, Y90 will be available five (5) days a week, Monday through Friday. Nuclitec has two manufacturing days: Wednesday and Friday.

The ordering process will remain the same. All Y90 orders will need to be placed with Zevalin Support Services at least one (1) week in advance of the Y90 treatment date.

Zevalin® Support Services phone number is 866-298-8433, Option 1.

Nuclitec Y90 is shipped and routed across US and its territories via FedEx under the usual Shipper's Declaration For Dangerous Goods.

Nuclitec will provide a Certificate of Analysis and the FedEx tracking number to each radiopharmacy in advance of the shipment date. Each shipment is expected to arrive at least one day **before** the calibration date. An example of the Certificate of Analysis is attached for your reference. It is batch specific.

All printed label data shown on the vial and package labels relating to product calibration and expiry will cite US EDT at 1200h.

Please note that the delivery vial is a 3 mL V-vial (Kimax) with substantially thicker walls. The vial is contained within a beta (& Bremstrahlung minimization) shield of clear acrylic with a removable top. The assembly is contained within a typical lead shield, as shown in the "thumbnail" picture below:






We look forward to working with Nuciltec and each radiopharmacy to provide flexibility in the treatment schedule for Zevalin® patients.

ZEVALIN Support Services PHONE: 866-298-8433 FAX: 877-264-8483 ZevalinSupport@sppirx.com





Nuclitec

Eckert & Ziegler Nuclitec GmbH

Gieselweg 1 38110 Braunschweig Deutschland www.nuclitec.de

Dr. Sven Dittrich Heat of Quality Control Medical

49 5307 932-261 49 5307 932-451 ven.darich@ezag.com

CERTIFICATE OF ANALYSIS

[90Y]Yttriumchloride in aqueous solution (0.04 m HCI)

Batch: YKD09

Characteristics	Requirement	Result
General data		1
Production date		2010-04-30
Reference date (as required by the customer)	Qate, 12h	as required by the customer
Expiry date (as required by the customer)	Reference date + 8 days	as required by the customer
Chemical or Physical Properties / Methods		S. C.
Appearance/visual check	Clear and colourless	passes test
Identity of Y-90/LSC	corresponds to	passes test
Gamma-emitting impurities (by the specified date)/gamma-spectrometry	x 10 ⁻⁴ Bq/Bq Y-90 (by the specified date)	< 3 x 10 ⁻⁵ Bq/Bq Y-90 07 May 2010, 12h CEST*
Sr-90 activity related to the Y-90 activity (by the specified date)/LSC	1 x 10 ⁻⁵ Bq Sr-90/Bq Y-90 (by the specified date)	< 2 x 10 ⁻⁷ Bq Sr-90/Bq Y-90 07 May 2010, 12h CEST*
Volume (as required by the customer)/ weighing (related to the customer requirements)	0.02 – 5.0 ml 90 – 110 %	passes test
Activity (given at reference date as required by the customer)/long ation chamber (related to the customer requirements)	0.1 – 300 GBq 90 – 110 %	passes test
DTPA-Binding/ LC	> 95 %	99.7 %
Radiochemical parity TLC	> 99 %	99.8 %
Other Tests		
Visual inspection of visible particles	free of particles	free of particles
Microbiological Specifications		
Autoclave parameters: examination of the autoclave log sterilising temperature [°C], sterilising time [min]	≥ +121.0 °C ≥ 20 min	+121.0 °C 20 min

* CEST = Central European Summer Time (UTC + 2h)

Eckert & Ziegler Nuclitec GmbH

i. A. Dr. Sven Dittrich **Head of Quality Control Medical**

Geschäftsführer Joseph Hathcock Dr. André Heß Wilhelm Wartenberg Sitz Braunschweig Amtsgericht Braunschweig Reg.-Nr. HRB 4647 USt-IdNr. DE193324337

Bankverbindung Commerzbank Braunschweig Bankleitzahl 270 400 80 Konto-Nr. 530 2591